



NEW

Superficial Artery Puncture Trainer

Multiple training scenarios, faithful reproduction of vessel behavior



Product No. MW66 | Part No. 11459-000 a set of 5

FEATURES

- | 01 Realistic tactile feedback as the dialysis needle penetrates the wall of the superficialized artery
- | 02 Wearable at any desired location on the limb using the Velcro strap
- | 03 Wearable design supports simulation training in a clinical setting

SKILLS

Palpation of puncture sites

Needle insertion into the superficialized artery

Observation of blood flashback

Notes: This trainer is not designed to simulate blood circulation using a hemodialysis machine.

[Recommended Devices]

Dialysis needle 16-18G

Simulates realistic vessel rolling during needle insertion.



DESCRIPTIONS

SPECIFICATIONS

Dimensions: W15×D11×H3.5 cm

Weight: 170 g (puncture unit)

Materials: soft resin/hard resin

SET INCLUDES

5	puncture unit	5	syringe 20mL
5	strap unit	10	simulated blood sticks (swab type)
5	bottle for simulated blood	5	plastic beaker
5	stand for bottle of simulated blood	1	storage case

RELATED PRODUCTS

M99 Arterial Puncture Training Wrist

MW82 Femoral Arterial Puncture Trainer



KYOTO KAGAKU



rw-kyoto@kyotokagaku.co.jp
www.kyotokagaku.com

June, 2026