MW76 11480-000



Automated Birthing Adapter "SAKUYA" for MW48

Consistent and Standardized Birth Assistance Training Via Mechanical Control







FEATURES

- 1 | Adjustable parameters for crowning-to-delivery time and a number of bearing-down efforts
- 2 | Standardized training made possible by mechanically controlled fetal descent and rotation
- 3 | Independent operation without assistants for limited-resource environments

SAKUYA is compatible with the following products

MW48 Perinatal Full-body Model "Hana" MW36 Delivery Assistance Simulator

SKILLS

Assistance with Delivery

- | Delivery Method: Supine Position
- | Observation of Crowning / Exposure
- | Protection of the Perineum / Anus
- | Fetal Position: Normal
- | Vacuum extraction
- | Delivery of the Baby
- | Delivery of the Placenta
- | Checking of the Amniotic Membranes

Neonatal Care

| Palpation of the Anterior and Posterior Fontanelles


Ligation and Cutting of the Umbilical Cord



F KEY FEATURES

Customizable fetal rotation direction, crowning repetitions, and total delivery time

The system features both automatic and manual modes, with simple operation via a touchscreen interface.



DESCRIPTIONS

SPECIFICATIONS

Size: W40x D36 x H36 cm W 15.6 x 14.2 D x H 14.2in. Weight: 3.6 kg / 7.9 lbs

1 power cable

lubricant (500ml)

MATERIALS

RECOMMENDED DEVICES

SET INCLUDES

- 1 abdominal unit
- 1 control box
 1 rotation unit
 - ition unit 1 instruction manual

1 fetal model

REPLACEMENT PARTS

OPTIONAL PARTS