




Automated Birthing Adapter "SAKUYA" for MW48

Consistent and Standardized Birth Assistance Training Via Mechanical Control




FEATURES

- 1 Adjustable parameters for crowning-to-delivery time and a number of bearing-down efforts
- 2 Standardized training made possible by mechanically controlled fetal descent and rotation
- 3 Independent operation without assistants for limited-resource environments

SAKUYA is compatible with the following products

MW48 Perinatal Full-body Model "Hana"
MW36 Delivery Assistance Simulator



SKILLS

Assistance with Delivery

- | Delivery Method: Supine Position
- | Observation of Crowning / Exposure
- | Protection of the Perineum / Anus
- | Fetal Position: Normal
- | Vacuum extraction
- | Delivery of the Baby
- | Delivery of the Placenta
- | Checking of the Amniotic Membranes


Neonatal Care

- | Palpation of the Anterior and Posterior Fontanelles
- | Ligation and Cutting of the Umbilical Cord

KEY FEATURES

Customizable fetal rotation direction, crowning repetitions, and total delivery time

The system features both automatic and manual modes, with simple operation via a touchscreen interface.



DESCRIPTIONS

SPECIFICATIONS

Size: W40x D36 x H36 cm
W 15.6 x 14.2 D x H 14.2in.

Weight: 3.6 kg / 7.9 lbs

SET INCLUDES

- | | |
|------------------|----------------------|
| 1 abdominal unit | 1 power cable |
| 1 control box | 1 lubricant (500ml) |
| 1 rotation unit | 1 instruction manual |
| 1 fetal model | |

MATERIALS

RECOMMENDED DEVICES

REPLACEMENT PARTS

OPTIONAL PARTS