


Smart Nursing and Emergency Care Torso JANA

Train with the Best
Empowering Caregivers of Medically Dependent Children



Patient handling

Airway suction: tracheostomy opening

Feeding tube management: oral / nasal / gastrostomy opening

Percutaneous Endoscopic Gastrostomy (PEG) care

Bag-Valve-Mask ventilation

Chest compressions and BLS skills assessment

Development of Educational Materials for Problem Resolution

1. This model is lightweight and highly portable, enabling practice among family members and multidisciplinary healthcare professionals.
2. Repetitive training can be conducted using this simulator [rather than on the child themselves].
3. Clear representation and visualization of the cannula site facilitates a better understanding of the principles underlying its management and care.
4. The model allows for practice of caregiving tasks associated with child-rearing, such as holding the child, replacing securing straps, and observing safety precautions during bathing.

Smart Nursing and Emergency Care Torso JANA


MW77 | 11481-000

■ Set includes

1	pediatric manikin
1	tracheostomy tube (sample)
1	gastrostomy tube (sample)
1	lubricant for catheterization
1	storage case
1	usage guide
1	instruction manual

■ Specifications

Size	W22 x D16 x H38 cm W8.7 x D6.3 x H15 in.
Weight	2.8 kg / 6.2 lbs
Materials	soft resin / hard resin
Application	"CPR add-on kit Student" Free download on App Store (iOS) Google Play Store (Android)



Designed to simulate the condition of a developing neck which lacks stability.

■ Recommended devices

tracheostomy tube : 18 Fr. (I.D. 4.0 mm O.D 6.0 mm)
gastrostomy tube : 14 Fr.

■ Replacement Parts



Lubricant
(for airway suction)
11229-050

■ Related products



MW21B
NCPR Simulator Plus II



MW71
Pediatric Physical
Examination Simulator

2025.05



<https://www.kyotokagaku.com/jp/>

rw-kyoto@kyotokagaku.co.jp

All American regions

Kyoto Kagaku America Inc. / Show Room

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA

Telephone: +1-310-325-8860 Fax: +1-310-325-8867

Email : info@kkamerica-inc.com

Europe, Russia & Africa

Kyoto Kagaku Europe GmbH / Training Room

De-Saint-Exupery-Str.10, 60549 Frankfurt, Germany

Telephone: +49-69-5060-28160

Asia and other regions

Kyoto Head Office / Factory / Show Room

15 Kitane-koya-cho Fushimi-ku Kyoto, Japan 612-8388

Telephone: +81-75-605-2510 Fax: +81-75-605-2519