

MW80



11484-000,-M,-D



NEW

Digital EAR Examination Simulator II

Digital system allows for expansion of case library of ear-drum images

Digital EAR Examination Simulator featuring image importer for educators to add a variety of different case images

Product Supervision

Japan Society for Medical Education
with the cooperation of Michihiko Sone
Professor- Nagoya university Graduatr
school of Medicine



FEATURES

- 1 | Custom fundus images can be programmed into the simulator via USB
- 2 | Ear canals come in 2 sizes - normal and stenosis
- 3 | Work with an actual otoscope
- 4 | The manikin neck can be tilted
- 5 | Training is possible with both left and right ears

CASES

normal tympanic membrane 1	acute suppurative otitis media (AOM)	otitis media with effusion
normal tympanic membrane 2	aerotitis media after myringotomy	retracted tympanic membrane
normal tympanic membrane 3	cholesteatoma 1	serous otitis media (SOM)
normal tympanic membrane 4	cholesteatoma 2	traumatic tympanic membrane perforation 1
normal tympanic membrane (pediatric)	cholesteatoma 3	traumatic tympanic membrane perforation 2
normal tympanic membrane (pediatric)	cholesteatoma 4	tuberculous otitis media
cancer of external auditory canal-adenoid cystic carcinoma	cholesterol granuloma of the middle ear	tympanic membrane atelectasis
cancer of external auditory canal-squamous cell carcinoma	chronic otitis media	tympanosclerosis
cerumen block	chronic otitis media with perforation 1	cochlear otosclerosis
cerumen	chronic otitis media with perforation 2	congenital ossicular anomaly - malleus bar
external auditory canal cholesteatoma	chronic otitis media with perforation 3	adhesive otitis media (pediatric)
external auditory canal exostosis	chronic otitis media with perforation 4	chronic otitis media with perforation (pediatric)
medial meatal fibrosis	eosinophilic otitis media	keratin cyst of the tympanic membrane (pediatric)
traumatic tympanic membrane perforation - post myringoplasty	glomus tympanicum tumor	otitis media with effusion (pediatric)
seborrheic keratosis of the external canal	large tympanic perforation	tympanic membrane atelectasis (pediatric)
acute otitis media	lateral tympanoplasty	postoperative ear (cholesteatoma)
	mucoid otitis media (MOM)	tympanostomy tube (pediatric)

KEY FEATURES



Selection by the button The back screen enables the instructor to monitor students and track the cases during examinations

SKILLS

- | Ear examination with a real otoscope
- | Ear wax and foreign body removal

50 Cases
Add +
Your Own Cases!

DESCRIPTIONS

SPECIFICATIONS

Size: 42 x 22 x 39 cm / 16.5 x 8.7 x 15.3 in.

Weight: 2.5kg / 5.5lbs

Power Consumption : 2.5W

SET INCLUDES

- 1 | manikin head and shoulder
- 1 | types of ear units -1 pair each
 - exam ear (normal)
 - exam ear (stenosis)
 - FB ear (normal)
 - FB ear (stenosis)

SET VARIATIONS

- 11484-000 MW80 light skin tone
- 11484-000-M MW80 medium skin tone
- 11484-000-D MW80 dark skin tone

REPLACEMENT PARTS

- 11391-010 ear units and foreign body set
 - 4 types of ear units - 1 pair each
 - ear exam (normal), ear exam (stenosis), FB ear (normal), FB (stenosis)
 - 1 set of foreign body removal kit
 - simulated ear wax, small objects and a spoon