

# MW83 12-Lead ECG Electrode Placement Trainer

An effective tool designed for fair skills assessment and examination under standardized conditions.

## Features

- **Trainer for sequential 12-lead ECG skills**, including intercostal space palpation and precise electrode placement.
- **Full-body manikin for developing non-technical skills**, including patient dignity and communication.
- **Anatomically accurate rib markings** to support skills assessment.
- **Modular design for easy transport and storage.**



\* The trainer has no electronic components and does not generate ECG signals.

## Skills

### Objective Skills Assessment for 12-Lead ECG

Assessment topics include, but are not limited to:

- Disinfection / Identification of appropriate electrode sites through intercostal space palpation
- Application of conductive paste / Electrode placement
- Consideration for patient comfort and dignity



# 12-Lead ECG Electrode Placement Trainer

MW83

11487-000

### Set Includes

manikin (upper body, legs)	1
chest skin (with rib marking)	1
pants	1
instruction manual	

\* ECG electrodes or monitor are not included in the set.

### Specifications

	full body adult manikin
dimensions	W55 x D18 x H160 cm
weight	8.2 kg
materials	soft resin



Replacement parts



Chest skin (with rib marking)  
11487-100

Optional parts



Chest skin (plain)  
11487-110

Related products



34 lung sound examples recorded from real patients


Lung Sound Auscultation Trainer "LSAT" ver.2  
MW28 11407-999



Covering essential skills for general physical assessment and vital signs measurement

Physical Assessment Simulator "Physiko Plus"  
M55B 11354-399/200

Intensive Pre-Exam Training with Physiko!



- ✓ Real-time visual feedback on correct electrode placement
- ✓ Dynamic ECG display on the monitor
- ✓ Twelve patient scenarios, including normal rhythm and myocardial infarction
- ✓ Connectable to standard ECG monitors



## Preparing OSCE?

### Looking for better tools for OSCE at your institution?

Our representatives are ready to listen to your needs and recommend the optimal combination of simulators. Contact us anytime.

**KYOTO KAGAKU**

<https://www.kyotokagaku.com/>  
[rw-kyoto@kyotokagaku.co.jp](mailto:rw-kyoto@kyotokagaku.co.jp)

2026.02

#### Kyoto Head Office / Factory / Show Room

15 Kitaneokoya-cho Fushimi-ku Kyoto, Japan 612-8388  
Telephone:+81-75-605-2510 Fax:+81-75-605-2519

#### Asia and other regions

##### Overseas Div. / MEM Training Center

Second fl. NREG Hongo-3-chome bldg., 3-26-6 Hongo,  
Bunkyo-ku, Tokyo 113-0033, Japan  
Telephone:+81-3-3817-8071 Fax:+81-3-3817-8075

#### All American regions

##### Kyoto Kagaku America Inc. / Show Room

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA  
Telephone:+1-310-325-8860 Fax:+1-310-325-8867  
Email : info@kkamerica-inc.com

#### Europe, Russia & Africa

##### Kyoto Kagaku Europe GmbH / Training Room

De-Saint-Exupery-Str.10, 60549 Frankfurt, Germany  
Telephone:+49-69-5060-28160