


# Clinical Assessment Kit “Cyanosis”, “Surgical Wound”

Putting on and sticking on  
Realistic patient setting in an instant

## Features

- | A handy versatility enhancer for existing simulators
- | Simple, tool-free and adhesive-free application
- | Peel-off colored “cyanosis” layer for flexible training scenarios
- | “Surgical Wound” kit featuring four reusable wound pad types (large/small, with/without sutures)



\*The above photos are to show the usage.

### Cyanosis

Peel-off colored layer

lips

Make appearance of cyanosis by applying the cyanosis painting on the lips and nails of simulators. It dries in 1-5 minutes using a dryer.  
A painting brush included

hand nails

toe nails

\*The photos here are preliminary images

### Surgical Wound Pad

Reusable

Durable silicone pad for quick, smooth attachment and removal  
A set of four wound types covering diverse training scenarios

Large: 13.5cm      Small: 8cm

## Skills

- Clinical interpretation and communication of findings
- Visual assessment of cyanosis
- Observation and assessment of surgical wounds

## Recommended Simulators to Apply

- MW25 Patient Care Simulator "Yaye"
- MW26 Basic Patient Care Simulator "Keiko"

\*The cyanosis kit is not compatible with MW22 Cherry\*

# Clinical Assessment Kit

MW85 | 11489-100 | Cyanosis  
 MW85 | 11489-000 | Surgical Wound Pad



**Set Includes** "Cyanosis"

cyanosis tube	4
painting brush	1
storage bag	4
instruction manual	

**Specifications**

	Cyanosis tube
contents	10g (per tube)



**Set Includes** "Surgical Wound Pad"

surgical wound pad (four types)	4 (1 each)
storage case	1
instruction manual	

**Specifications**

	Surgical wound pad
dimensions	[large pad] W13.5×D0.3×H2.5cm [small pad] W 8×D0.3×H2cm

Related products

## Pursuing Greater Realism



simulated pus  
(set of 2, 100 ml each)  
11446-020



simulated exudate  
(set of 2, 100 ml each)  
11446-030



simulated blood  
(set of 2, 100 ml each)  
11446-040

When applied to the surgical pad, these simulated wound fluids create a true-to-life appearance of hemorrhage or exudate.

## Use the MW85 Clinical Assessment Kit on:



An integrated patient care simulator for a wide range of trainees, from undergraduates to professionals.

MW25 Patient Care Simulator "Yaye"  
11404-000



A basic patient care simulator for essential nursing and caregiving skills.

MW26 Basic Patient Care Simulator 1405-000

\* The cyanosis kit is not compatible with MW22 Cherry \*



<https://www.kyotokagaku.com/jp/>  
[rw-kyoto@kyotokagaku.co.jp](mailto:rw-kyoto@kyotokagaku.co.jp)

2026.03 KK

Kyoto Head Office / Factory / Show Room  
 15 Kitane-koya-cho Fushimi-ku Kyoto, Japan 612-8388  
 Telephone:+81-75-605-2510 Fax:+81-75-605-2519

Asia and other regions  
 Overseas Div. / MEM Training Center  
 Second fl. NREG Hongo-3-chome bldg., 3-26-6 Hongo,  
 Bunkyo-ku, Tokyo 113-0033,  
 Telephone:+81-3-3817-8071 Fax:+81-3-3817-8075

All American regions  
 Kyoto Kagaku America Inc. / Show Room  
 3109 Lomita Boulevard, Torrance, CA 90505-5108, USA  
 Telephone:+1-310-325-8860 Fax:+1-310-325-8867  
 Email : info@kkamerica-inc.com

Europe, Russia & Africa  
 Kyoto Kagaku Europe GmbH / Training Room  
 De-Saint-Exupery-Str.10, 60549 Frankfurt, Germany  
 Telephone:+49-69-5060-28160