

NEW

Newborn Bathing and Nursery Care Model II Ken and Rica

Adorably lifelike, expressive newborns for a wide range of training scenarios.



MW87
MW88

11491-000 Ken (male)
11492-000 Rica (female)

FEATURES

- | 01 Flexible, non-rigid neck
- | 02 A one-piece, seamless waterproof body with soft earlobes supports baby-bathing training
- | 03 Palpable anterior and posterior fontanelles

SKILLS

- Bathing
- Process of physical assessment
- Assessment of a newborn's overall condition
- Umbilical cord care
- Palpation of fontanelles
- Assessment of head control (neck stability)

The mouth can accommodate both breastfeeding and bottle-feeding nipples.



CASES/PATHOLOGIES/ORGANS

Anterior fontanelle
Posterior fontanelle
Umbilical cord

DESCRIPTIONS

Specifications

Full-body newborn manikin (male /female)
Dimensions : 50cm tall
Head circumference : 34cm

Weight : 約3.6kg
Materials : Soft resin

Set Includes

1 newborn manikin (male or female)
1 umbilical cord



KYOTO KAGAKU



rw-kyoto@kyotokagaku.co.jp