


MW8/MW8A Tube Feeding Chest Sheet

Instruction Manual

This sheet helps trainees anatomical understanding before the training of tube feeding.



● Cover the model with this sheet

Remove the skin before setting the sheet on the model.

Do not set the sheet directly on the skin so that the printed ink of the sheet marks on the skin and is irremovable.

● How to fix the sheet

Set the sheet with right and left button of the model.



Caution

● Handle the sheet with care.

Don't remove and handle the sheet forcibly. Handle the sheet with care to prevent from being damaged.

● This sheet may be cleaned with mildly soapy water or diluted detergent.

Never use ethanol or organic solvent like paint thinner so that it may cause the transformation and removing the print.

● Storage in a dark, cool space.

Do not fold the sheet and store it in the attached bag.

• For inquiries and service, please contact your distributor or KYOTO KAGAKU CO., LTD.

2026.08.07

 **KYOTO KAGAKU CO., LTD**



<http://www.kyotokagaku.com>



rw-kyoto@kyotokagaku.co.jp

Contact Address:

LA Office (for USA, CANADA and Mexico customers)

Kyoto Kagaku USA 3109 Lomita Blvd, Torrance, CA 90505

Phone: 310-325-8860 Fax: 310-325-8867

Main Office and Factory

15 Kitane-koya-cho Fushimi-ku Kyoto 612-8388, Japan

Phone : 81-75-605-2510

Fax : 81-75-605-2519