

PH-17 | 41317-000

## Water Body Phantom WAC

X-ray



Chest and abdomen phantom to help measure surface dose



SHOW MORE!



## FEATURES

- Water body phantom represents human chest and abdomen to serve as radiation absorber and scatterer.
- The phantom is a double walled cylindrical container with track shaped cross section.
- For chest, fill water between two walls.
- For abdomen, fill the both spaces.

## DESCRIPTIONS

## SET INCLUDES

- 1 body phantom
- 1 storage case

## SPECIFICATIONS

Phantom size:  
W30×D20×H45 cm  
W11.8×D7.9×H17.7 in


PH-14 | 41430-000-

## Acrylic Phantom XAC

X-ray



Slab phantoms for radiation absorption and scattering measurement



SHOW MORE!



## VARIATIONS

|        |           |                                        |
|--------|-----------|----------------------------------------|
| XAC-01 | 41430-000 | 30 × 30 × 0.1 cm/11.8 × 11.8 × 0.04 in |
| XAC-02 | 41431-000 | 30 × 30 × 0.2 cm/11.8 × 11.8 × 0.08 in |
| XAC-03 | 41432-000 | 30 × 30 × 0.3 cm/11.8 × 11.8 × 0.12 in |
| XAC-04 | 41433-000 | 30 × 30 × 0.4 cm/11.8 × 11.8 × 0.16 in |
| XAC-05 | 41434-000 | 30 × 30 × 0.5 cm/11.8 × 11.8 × 0.2 in  |
| XAC-08 | 41435-000 | 30 × 30 × 0.8 cm/11.8 × 11.8 × 0.3 in  |
| XAC-1  | 41436-000 | 30 × 30 × 1 cm/11.8 × 11.8 × 0.4 in    |
| XAC-2  | 41437-000 | 30 × 30 × 2 cm/11.8 × 11.8 × 0.8 in    |
| XAC-3  | 41438-000 | 30 × 30 × 3 cm/11.8 × 11.8 × 1.2 in    |
| XAC-4  | 41439-000 | 30 × 30 × 4 cm/11.8 × 11.8 × 1.6 in    |
| XAC-5  | 41440-000 | 30 × 30 × 5 cm/11.8 × 11.8 × 2.0 in    |
| XAC-8  | 41441-000 | 30 × 30 × 8 cm/11.8 × 11.8 × 3.1 in    |
| XAC-10 | 41442-000 | 30 × 30 × 10 cm/11.8 × 11.8 × 3.9 in   |

PH-16 | 41318-000,010 / 41319-000,010

## Contrast Detail Phantom

X-ray



Concentration resolution evaluation from two directions in plain X-ray



SHOW MORE!



## FEATURES

- Four types of phantoms with different sizes and target types
- Hole 15: 41318-000      - Hole 10: 41318-010
- Rod 15: 41319-000      - Rod 10: 41319-010

## DESCRIPTIONS

## SIZE

- Hole 15: 15 × 15 holes of depth range from 1.0 to 8.0 mm (0.4 to 3.1 in)
- Rod 15: 15 × 15 rods of height range from 1.0 to 8.0 mm (0.4 to 3.1 in)
- Hole 10: 10 × 10 holes of depth range from 1.0 to 5.5 mm (0.4 to 2.2 in)
- Rod 10: 10 × 10 rods of height range from 1.0 to 5.5 mm (0.4 to 2.2 in)

## SET INCLUDES

- 1 chart phantom
- storage case