

PH-17 | 41317-000

Water Body Phantom WAC

X-ray



Chest and abdomen phantom to help measure surface dose



SHOW MORE!



FEATURES

- Water body phantom represents human chest and abdomen to serve as radiation absorber and scatterer.
- The phantom is a double walled cylindrical container with track shaped cross section.
- For chest, fill water between two walls.
- For abdomen, fill the both spaces.

DESCRIPTIONS

SET INCLUDES

- 1 body phantom
- 1 storage case

SPECIFICATIONS

Phantom size:
W30×D20×H45 cm
W11.8×D7.9×H17.7 in


PH-14 | 41430-000-

Acrylic Phantom XAC

X-ray



Slab phantoms for radiation absorption and scattering measurement



SHOW MORE!



VARIATIONS

XAC-01	41430-000	30 × 30 × 0.1 cm/11.8 × 11.8 × 0.04 in
XAC-02	41431-000	30 × 30 × 0.2 cm/11.8 × 11.8 × 0.08 in
XAC-03	41432-000	30 × 30 × 0.3 cm/11.8 × 11.8 × 0.12 in
XAC-04	41433-000	30 × 30 × 0.4 cm/11.8 × 11.8 × 0.16 in
XAC-05	41434-000	30 × 30 × 0.5 cm/11.8 × 11.8 × 0.2 in
XAC-08	41435-000	30 × 30 × 0.8 cm/11.8 × 11.8 × 0.3 in
XAC-1	41436-000	30 × 30 × 1 cm/11.8 × 11.8 × 0.4 in
XAC-2	41437-000	30 × 30 × 2 cm/11.8 × 11.8 × 0.8 in
XAC-3	41438-000	30 × 30 × 3 cm/11.8 × 11.8 × 1.2 in
XAC-4	41439-000	30 × 30 × 4 cm/11.8 × 11.8 × 1.6 in
XAC-5	41440-000	30 × 30 × 5 cm/11.8 × 11.8 × 2.0 in
XAC-8	41441-000	30 × 30 × 8 cm/11.8 × 11.8 × 3.1 in
XAC-10	41442-000	30 × 30 × 10 cm/11.8 × 11.8 × 3.9 in

PH-16 | 41318-000,010 / 41319-000,010

Contrast Detail Phantom

X-ray



Concentration resolution evaluation from two directions in plain X-ray



SHOW MORE!



FEATURES

- Four types of phantoms with different sizes and target types
- Hole 15: 41318-000 - Hole 10: 41318-010
- Rod 15: 41319-000 - Rod 10: 41319-010

DESCRIPTIONS

SIZE

- Hole 15: 15 × 15 holes of depth range from 1.0 to 8.0 mm (0.4 to 3.1 in)
- Rod 15: 15 × 15 rods of height range from 1.0 to 8.0 mm (0.4 to 3.1 in)
- Hole 10: 10 × 10 holes of depth range from 1.0 to 5.5 mm (0.4 to 2.2 in)
- Rod 10: 10 × 10 rods of height range from 1.0 to 5.5 mm (0.4 to 2.2 in)

SET INCLUDES

- 1 chart phantom
- storage case