

NEW

Breast plate for Chest Phantom N-1 LUNGMAN

41337-090



Breast plate for Chest Phantom N-1 "LUNGMAN". This plate simulates the absorption of breast in female chest X-ray, and the influence of nipple on image quality.

Specification
Size : 47×11×55 cm
Weight : 5.2kg
Material : urethane based resin (specific gravity: 1.06)

Related Products

Pneumonia module for Chest Phantom N-1 LUNGMAN



Chest phantom close to the actual image of pneumonia.

Specification
Size : 13×16×19 cm
Weight : 2.2kg
Material : urethane based resin (specific gravity: 0.92)