

PH-79 | 41945-000

## X-Ray Training Phantom PBU-POSE


X-ray

ANTHROPOMORPHIC PHANTOMS

**For patient-friendly and accurate positioning learning**  
**Supports scenario based trainings including communication skills**



\*The wheelchair is not included.



SHOW MORE!



## FEATURES

- | Light weight and close-to-human articulation as well as anatomical landmarks for positioning and patient handling.
- | Soft touch of phantom's skin facilitates realistic simulation training scenario including communication skills.
- | Radiography with a lower irradiation reducing risk for trainees and stress on the device.
- | Enables training free from privacy concerns and inconveniences associated with use of standardized patients.

## APPLICATIONS

- | Patient positioning
- | Patient transportation
- | Plain radiography

## ANATOMY

## Skeletal system


skull, cervical spine, vertebrae, clavicles, scapulae, sternum, pelvis, lungs (no vessels), heart, kidneys, upper and lower arm bones, carpal, metacarpal, femur, kneecaps, lower leg bones, tarsi, metatarsals, phalanges

## Internal organs

trachea (up to 1st bifurcation), lungs (diaphragm only), heart, kidneys

## Land marks

external acoustic foramen, mastoid, seventh cervical vertebra, manubrium, xiphisternum, styloid process of radius, superior margin of the symphysis pubis, medial epicondyle of femur/epicondylus lateralis, patella, malleolus (internal condyle /external condyle), subcostal area, landmark on the body surface, trochanter, processus styloideus, ulnae



## SET INCLUDES

1 whole body phantom	1 set of assembly tools
1 head supporter	1 set of sample X-ray data (DVD)
1 pajamas	1 manual

## SPECIFICATIONS

## Phantom dimensions

height:	chest girth : 85cm	(thickness : 20cm)
165 cm/64.9in	33.5in	(thickness : 7.9cm)
weight:	waist girth : 75cm	(thickness : 19cm)
18 kg / 39.69 lb	29.5 in	(thickness : 7.5cm)

## MATERIALS

Soft fabrics : polyurethane foam (density 0.2)  
 Skeleton : epoxy resin (density 1.31)  
 Skull : urethane resin (density 1.12)

## Packing size:

W113 x D59 x H37 cm  
 W 44.5x D23.2 x H 14.6 in

