

US-25

41957-000




NEW

# Thyroid Biopsy Phantom

## Enhance confidence in performing safe and accurate thyroid biopsies

Visualization of the thyroid and targets in reference to sonographic landmarks.




### FEATURES


- 1 | Transparent and opaque biopsy blocks for introductory and advanced learners, featuring step-by-step acquisition of the biopsy skills.
- 2 | Targets in two sizes containing liquid for aspiration are embedded in the thyroid.
- 3 | Targets and carotids are colored to confirm successful sample collections as well as complications.

### SKILLS

- | Ultrasound examination on the neck
- | Ultrasound-guided FNA (Fine Needle Aspiration)



Internal structure



Target and the FNA needle

### SPECIFICATIONS

Size: W26 x D14 x H 12 cm  
W10.2 x D5.5 x H4.7 in.

Weight: 1.5kg / 3.3lbs

Target size: 5mm dia., 7mm dia.  
0.27in dia., 0.19in dia.

### SET INCLUDES

- 1 face and neck manikin
- 1 transparent biopsy block
- 1 opaque biopsy block
- 1 storage case
- 1 instruction manual

### DESCRIPTIONS

#### MATERIALS

Soft resin, hydrogel, hard resin  
Latex free

#### REPLACEMENT PARTS

41957-010 2 biopsy blocks  
(duo set :  
transparent + opaque each one)

Ultrasound scanner and needles are not included.