

US-10 US-10A US-10B US-10C US-10D



Female Pelvic Ultrasound Phantom

Practice both transvaginal and transabdominal ultrasound procedures in one

Newly added 10-week- pregnancy unit and IUD unit supports thorough training in Obs/gyn ultrasound

Product Supervision

Charlotte Henningsen, MS, RT(R), RDMS, RVT,
FSDMS, FAIUM
Chair & Professor - Sonography Department
Adventist University of Health Sciences



FEATURES

- 1 | Compatible with actual ultrasound transducer
- 2 | Five types of interchangeable phantom inserts with different pathologies
- 3 | Realistic view of female external genital
- 4 | Excellent ultrasound image quality

SKILLS













- | Handling and manipulation of transvaginal and transabdominal transducers
- | Interpretation of sonogram
- | Visualization and localization of anatomies and pathologies

DESCRIPTIONS

SET INCLUDES

	US-10A/41909-100	US-10C/41909-110	US-10/41909-000	US-10B/41909-200	US-10D/41909-210
	Female Pelvic Ultrasound Phantom (3units/7w)	Female Pelvic Ultrasound Phantom (3units/10w)	Female Pelvic Ultrasound phantom (2 units)	Obstetric Pelvic Ultrasound Phantom (7w)	Obstetric Pelvic Ultrasound Phantom (10w)
lower torso manikin	✓	✓	✓	✓	✓
ultrasound pathology unit	✓	✓	✓	—	—
ultrasound 7-week pregnancy unit	✓	—	—	✓	—
ultrasound 10-week pregnancy unit	—	✓	—	—	✓
ultrasound ectopic pregnancy unit	✓	✓	✓	✓	✓
IUD unit	—	—	—	—	—
storage case	2	2	1	1	1
talcum powder	✓	✓	✓	✓	✓
syringe	✓	✓	✓	✓	✓
instruction manual	✓	✓	✓	✓	✓

 KEY FEATURES

Five variations of training sets	Anatomy, findings and pathology	US-10A	US-10C	US-10	US-10B	US-10D
Normal pregnancy (7-week)	  [Longitudinal image]	✓	—	—	✓	—
Normal pregnancy (10-week)	  [Longitudinal image]	—	✓	—	—	✓
Ectopic pregnancy unit	 	✓	✓	✓	✓	✓
Pathological unit	  Anatomies Pathologies	✓	✓	✓	—	—
IUD unit	 	—	—	—	—	—
NEW Follicular phase unit	 	—	—	—	—	—

SPECIFICATIONS

Torso size: W34 x D33 x H24 cm
W13.5 x D13 x H9.5 in.

weight: 7.2kg / 15.9 lbs

DESCRIPTIONS**MATERIALS**

Soft resin, Polyurethane elastomer
Latex free

* Ultrasound scanner is not included.

REPLACEMENT PARTS

41909-030 ultrasound pathology unit
41909-040 ultrasound ectopic pregnancy unit
41909-050 ultrasound first trimester pregnancy unit (7W)
41909-060 skin for US-10 / 10A / 10B
41909-070 ultrasound first trimester pregnancy unit (10W)
41909-080 IUD unit
41909-090 follicular phase unit