US-11

Scrotal Ultrasound Phantom

Master techniques for a safe and painless examination

Hands-on training to accurately diagnose scrotum abnormalities

Product Supervision

Charlotte Henningsen, MS, RT(R), RDMS, RVT, ESDMS, FAIUM Chair & Professor - Sonography Department Adventist University of Health Sciences

FEATURES

- 1 | Two types of interchangeable scrotums: normal and with pathologies
- 2 | Excellent ultrasound image quality


SKILLS

Patient positioning and preparation for examination Scrotal ultrasound screening | Visualization of testicular cancer

KEY FEATURES


Normal phantom


- scrotum, testicle, epididymis and penis



Pathological phantom

-scrotum, testicle, epididymis and penis -testicular cancer (each one in left and right testis, 10mm dia.)





SPECIFICATIONS

Size: 34 x 33 x 24 cm/13.4 x 13 x 9.4 inch SET INCLUDES

Weight: 3.6kg / 7.9lbs

DESCRIPTIONS

REPLACEMENT PARTS

41911-010 1 normal scrotum unit 41911-020 1 pathological scrotum unit

1 pelvic body

- normal scrotum 1
- 1 pathological scrotum 1
- storage case
- 1 instruction manual talcum powder

KYOTO KAGAKU