

US-1



Ultrasound Examination Training Phantom "ECHOZY"

Abdominal ultrasound phantom without pathologies

Unique high-fidelity abdominal ultrasound phantom that facilitates effective training compatible with your own clinical devices.

Product Supervision

Dr. Hitoshi Asai, Director
Osaka Kyoiku University, Health Administration Center
Dr. Shigeru Nakamura
Nagayoshi General Hospital, Clinical Examination
Department



Full set includes
anatomical model
"ECHO-ZOU"

FEATURES


- 1 | Detailed hepatobiliary, pancreatic and other abdominal anatomy
- 2 | Eight Couinaud's hepatic segments can be localized.

SKILLS

- | Basics of abdominal sonography:
 - Cross sections and sonographic anatomy
 - Sonographic demonstration of each individual organ
 - Localization of Couinaud's hepatic segments

Anatomy

- | Liver
(segmental anatomy, portal and hepatic venous systems, ligamentum teres and ligamentum venosum)
- | Biliary tract
(gallbladder, cystic duct, intrahepatic and extrahepatic bile ducts)
- | Pancreas (pancreatic duct)
- | Spleen / kidneys
- | Detailed vascular structures
(aorta, vena cava, celiac artery and its branches, portal vein and its branches, superior mesenteric vessels, renal vessels, and more)



DESCRIPTIONS

SPECIFICATIONS

Size: 29 x 19 x 31 cm/11.4 x 7 x 12.2 inch Weight: 12kg / 26.5lbs

MATERIALS

Polyurethane elastomer
Latex free

SET INCLUDES

US-1 (41900-010) "ECHOZY"	US-1 (41900-000) "ECHOZY" full set
1 ultrasound phantom "ECHOZY"	1 ultrasound phantom "ECHOZY"
1 set positioning pillow	1 anatomical model "ECHO-ZOU"
1 talcum powder	1 positioning pillow
1 storage case	1 tutorial manual (DVD)
	1 talcum powder
	2 storage cases