



Ultrasound Examination Training Phantom ECHOZY

Abdominal ultrasound phantom without pathologies

Unique high-fidelity abdominal ultrasound phantom that facilitates effective training compatible with your own clinical devices

Product Supervision

Dr. Hitoshi Asai, Director
Osaka Kyoiku University, Health Administration Center
Dr. Shigeru Nakamura
Nagayoshi General Hospital, Clinical Examination Department



Full set includes anatomical model ECHO-ZOU (11224-000)

FEATURES



- 1 | Detailed hepatobiliary, pancreatic and other abdominal anatomy
- 2 | Couinaud's eight hepatic segments can be localized

SKILLS

- | Basics of abdominal sonography:
 - Cross sections and sonographic anatomy
 - Sonographic demonstration of each individual organ
 - Localization of Couinaud's hepatic segments

Anatomy

- | Liver (segmental anatomy, portal and hepatic venous systems, ligamentum teres and ligamentum venosum)
- | Biliary tract (gallbladder, cystic duct, intrahepatic and extrahepatic bile ducts)
- | Pancreas (pancreatic duct)
- | Spleen / kidneys
- | Detailed vascular structures (aorta, vena cava, celiac artery and its branches, portal vein and its branches, superior mesenteric vessels, renal vessels, and more)



DESCRIPTIONS

SPECIFICATIONS

Size: W29 x D19 x H31 cm / W11.4 x D7 x H12.2 in. Weight: 12kg / 26.5lbs

MATERIALS

Polyurethane elastomer
Latex free

SET INCLUDES

US-1 (41900-010) ECHOZY	US-1 (41900-000) ECHOZY full set
1 ultrasound phantom ECHOZY	1 ultrasound phantom ECHOZY
1 set positioning pillow	1 anatomical model ECHO-ZOU
1 tutorial manual (DVD)	1 positioning pillow
1 talcum powder	1 tutorial manual (DVD)
1 storage case	1 talcum powder
	2 storage cases