



15 Kitanechoya-cho Fushimi-ku

Kyoto, Japan 612-8388

Telephone : +81-75-605-2510

Fax : +81-75-605-2519

Dear customer,

Note

Thank you for choosing our ultra sound phantom SPACEFAN-ST.

There may be a case that the surface of SPACEFAN-ST is in a wet condition. The ultra sound image will not be affected. In such case, wipe the surface with medical use alcohol. **Do not wipe the surface with strong solvent such as benzene.**

If you have any questions regarding this matter, please do not hesitate to contact us at any time.

Faithfully yours,

Kyoto Kagaku Co., Ltd.

Email: rw-kyoto@kyotokagaku.co.jp