




# Pediatric FAST/Acute Abdomen Phantom

## Ultrasound training for pediatric FAST procedures

This world's first pediatric ultrasound torso phantom provides opportunities of hands-on training in ultrasound examination of a two-year-old patient

### Product Supervision

Shunsuke Nosaka, MD, PhD  
 Director of Radiology  
 National Medical Center for Children and Mothers



### FEATURES

- 1 | Thoracoabdominal organs, bone structure and free-fluid for practicing FAST examination procedures
- 2 | Commonly seen pathologies in pediatrics
- 3 | Acquire fundamental pediatric abdominal ultrasound skills, as well as FAST/Acute sonography skills


### SKILLS

- | FAST procedures:  
Internal hemorrhage at perihepatic, perisplenic, pelvis, and pericardium area
- | Sonography for acute pediatric patient:  
Bowel intussusception, appendicitis and biliary dilatation

### KEY FEATURES



Pleural hemorrhage



Hydronephrosis



Appendicitis

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W41 x D17 x H15 cm      Weight: 7kg / 15.4lbs  
 W16.1 x D6.7 x H5.9 inch

#### SET INCLUDES

- 1 | ultrasound phantom
- 1 | talcum powder
- 1 | storage case
- 1 | tutorial manual (DVD)
- 1 | instruction manual

#### MATERIALS

Polyurethane elastomer  
 Latex free

#### REPLACEMENT PARTS

#### RECOMMENDED DEVICES

\* Ultrasound scanner is not included.