

MW26




# Basic Patient Care Simulator "Keiko"

## Meet Keiko :

**Patient simulator of geriatric care for trainees including home-care nurses and caregivers**

Keiko simulates patients under long-term cares in a variety of health care settings



**Feeding**

- | Tube feeding (oral/nasal)

**Exercise**

- | Body positioning
- | Patient transfer

**Wound Care**

- | Pressure sore care

**Pain Management**

- | Hot and Cold Compressions

**Elimination**

- | Enema
- | Urinary catheterization
- | Stoma care
- | Bedpan practice

**Hygiene**

- | Oral care
- | Bed bathing
- | Change of clothes
- | Perineal care
- | Hair care\*
- | Grooming and appearance

**Respiratory Care**

- | Oxygen inhalation
- | Airway suction
- | Postural drainage
- | Tracheotomy care

**Others**

- | Postmortem care

**Drug Administration**

- | Suppository

\* Optional



**KEY FEATURES**

**Basic Patient Care**

**Oral Care**

- Oral care on the bed and wheel chair
- Special care for patients with feeding tubes



- | Oral assessment
- | Teeth brushing
- | Denture care

**Pressure sores care**

- Learn stages of pressure sores
- Pressure sore care



- | Stage 2: Scapula
- | Stage 3: Great trochanter
- | Stage 4: Sacrum bone
- | Unstageable: Heel


**Perineal care**

- Anatomically correct genitals for detailed care



**Body positioning and transfer**

- Communication
- Correct manual handling techniques
- Patient safety



- | Changing body positions
- | Supporting patient transfer

**Medical Understanding**

**Airway Suction**

- Understanding the procedure



- | Airway suction: oral, nasal and tracheotomy (\*tube insertion only)
- | Tube feeding: oral, nasal and PEG (\*tube insertion only)

**Daily care with medical understandings**



- | Stoma care
- | Bed bathing for the patients with central venous port
- | Perineal care for the patients with urinal catheters



Photo: Kyoto College of Nursing

\* A wheelchair is not included.

**DESCRIPTIONS**

**SPECIFICATIONS**

Size: 149 cm/58.7 inch      Weight: 15kg/ 33lbs

**SET INCLUDES**

- 1 manikin
- 1 denture
- 1 wig
- 1 oral cover for suctioning
- 1 female genitalia

**MATERIALS**

Latex free  
Soft resin

- 1 male genitalia
- 1 under pants
- 1 wrench
- 1 instruction manual

**REPLACEMENT PARTS**

- 11405-010 female genitalia
- 11405-020 male genitalia
- 11401-080 denture
- 11229-060 oral cover for suctioning
- 11229-050 lubricant for training models
- 11381-800 lubricant for catheterization

**OPTIONAL PARTS**

- 12746-020 tracheotomy tube
- 11251-320 simulated suppository set
- 11251-330 simulated feces set
- 11405-030 pajama
- 11404-030 washable wig (elderly)
- 11229-070 sputum

**RECOMMENDED DEVICES**

- Nasogastric tube feeding:
  - Nasal catheter 12fr -14fr
  - Gastrostomy catheter for replacements 20fr
- Indwelling Urinary Catheter:
  - 16fr (male)/ 14fr (female)