



Obstetric Training Simulator -complete set-

Comprehensive maternal body set for skills training in vaginal examinations during labor and delivery



FEATURES

1 | The external genitalia is made of innovative material that excels in elasticity durability resilience for effective and repetitive training.

SKILLS

- | Cervical examination
- | Delivery assistance
- | Perineorrhaphy (perineal repair)

KEY FEATURES


The external genitalia is made of innovative material that excels in elasticity durability resilience for effective and repetitive training.





LANDMARKS

velamen (amniotic sac) vaginal canal ischial spine pubic symphysis соссух anterior fontanelle





5 types of cervical dilation models and fetal descent positions can be interchanged with a single action.

- | Cervical scoring in the first stage of labor
 - -Determination of station using landmarks of the ischial spine and the pubic joint to assess the descent of the body
 - -Assessment of cervical dilation
- Speculum examination
- Observation and cervical cytology
- Urinary catheterization









- Delivery assistance for different delivery positions (dorsal, all-fours)
- | Perineal protection
- | Delivery of fetus

(forceps delivery, vacuum extraction)

- | Clamping, tying and cutting of the umbilical cord
- Delivery of the placenta
- | Insertion of urinary catheter













Repetitive suture training for the first-degree perineal tears.

First-degree perineal tear repair

- | Single interrupted suture
- | Vertical mattress suture



DESCRIPTIONS

SPECIFICATIONS

Size: W74 x D44 x H33.5cm/ W29.1x D17.3 x H13.2 inch

SET INCLUDES

- 1 maternal body
- genitalia unit (for pelvic examination and delivery assistance)
- talcum powder
- 3 lubricant (150ml)
- instruction manual

Weight: 11.2kg/ lbs

Parts for pelvic examination:

- 1 vaginal unit set of cervical dilation inserts (5 variations)

 - 1 holder for the cervical dilation inserts

MATERIALS

Soft resin

Parts for Delivery Assistance:

- 1 fetus
- 3 umbilical cords for omphalotomy
- trav
- 5 velamen sheets



