

MW43

NEW

Abdominal Simulator for Examination

True-to-life representation of tactile features of abdominal organs, which are convincing enough to serve as a benchmark in cultivating the sense of touch.

Product Supervision

Takahiro Amano, M.D.,
Senior Advisor for Academic Affairs and Professor
International University of Health and Welfare,
Honorary Director Sannno Medical Center, Japan


Kyoko Arai
RN, PHN, PhD, Professor,
School of Nursing, University of Shizuoka, Japan



Interchangeable organs



Supple, resilient and delicate texture



Respiratory motion

FEATURES

- 1 | Supple, resilient and delicate textures which allow both shallow and deep palpation
- 2 | Interchangeable organs
- 3 | Respiratory motion
- 4 | Abdominal auscultation

SKILLS

- | Visual Inspection
- | Auscultation
 - renal artery
 - abdominal aorta
 - iliac artery
 - bowel sounds
- | Palpation
 - liver, spleen, kidney, uterus and lymphoma
- | Percussion
 - liver, spleen and rib

KEY FEATURES

Bony landmarks


- | Pelvis
- | Ribs
- | Lower rib cage
- | Costal margin
- | Xiphisternum
- | Pubic crest
- | Anterior superior iliac spines

Liver

- | Normal
- | Chronic hepatitis
- | Early late cirrhosis
- | Cirrhosis

Inguinal area

- | Normal
- | Lymphoma



Spleen

- | Normal
- | Enlarged
- | Slightly enlarged

Kidney

- | Normal
- | Hydronephrosis
- | Polycystic kidney

Uterus

- | Normal
- | Myoma

Bruit and sounds

- Renal artery
- | Abdominal aorta
- | Iliac artery
- | Bowel sounds

DESCRIPTIONS

SPECIFICATIONS

Size: 40x80x20cm/15.7x31.5x7.9 inch

Weight: 14kg/ 30.9lbs

MATERIALS

Soft resin
Latex free

SET INCLUDES

- | | |
|---|---|
| 1 male torso unit | 2 inguinal lymph nodes (normal, lymphoma) |
| 1 talcum powder | 2 uterus (normal, fibroid) |
| 1 set of the respiratory motion control unit | 1 simulated stethoscope |
| 2 kidneys (hydronephrosis, cysts) | 1 controller |
| 2 spleens (slightly enlarged, enlarged) | 1 AC adapter |
| 4 livers (precirrhosis, cirrhosis, chronic hepatitis, normal) | 4 storage cases |
| | 1 instruction manual |