




FAST/Acute Abdomen Phantom FAST/ER FAN

Conduct "FAST" assessment to save a trauma patient!

FAST/ER FAN provides simulated training in FAST (Focused Assessment with Sonography for Trauma); an ultrasound examination directed at detection of the presence of free fluid in trauma patients

Product Supervision

Junji Machi, MD, PhD
University of Hawaii at Manoa and Kuakini Medical Center



FEATURES

- 1 | An innovative phantom for repetitive training of FAST as an adjunct to the ATLS primary survey
- 2 | Pathologies include cholecystitis, an aortic aneurysm and a lesion on the colon




SKILLS

- | Internal hemorrhage at perihepatic, perisplenic, pelvis and pericardium area
- | Sonography for acute patients

FAST



Acute abdomen



DESCRIPTIONS

SPECIFICATIONS

Size: W62 x D30 x H25 cm Weight: 30kg / 66lbs
W24.4 x D11.8 x H9.8 in.

SET INCLUDES

- 1 | ultrasound phantom with carrying base
- 1 | tutorial manual (DVD)
- 1 | talcum powder
- 1 | instruction manual

MATERIALS

Polyurethane elastomer
Latex free

REPAIR

41903-000-99 recoating of US-5

*Ultrasound scanner is not included.