Caution

Do not mark on the product with a pen or leave printed materials in contact with its surface. Ink marks may be absorbed and become irremovable.


PH-80

Daily QA Phantom WEM Aqua Slab

Instruction Manual

Contents

- Set Includes, General Precaution------P.1
- Reference--------------P.2



KYOTO KAGAKU CO.,LTD

Before You Start

Set Includes/ Specifications



Material
Polyurethane (Aqua Slab)

41948-000 without attachment plate
Dimensions W20 × D20 × H20 cm
Weight 6.4 kg
41948-100 with attachment plate
Dimensions W20 × D20 × H21.5 cm
Weight 6.5 kg

Before first use, verify that all components listed below are included.

41948-000 without attachment plate
Phantom 1
41948-100 with attachment plate
Phantom with an acrylic plate 1
Attachment plate 1


General Precautions

●Handle with Care Handle the product and its components with care. Dropping or applying strong impact may cause damage.	●Clean the phantom with a soft cloth moistened with water or mild detergent. Never wipe the phantom or components with thinner or other organic solvents.
●Store the phantom at room temperature, away from heat, moisture, and direct sunlight to avoid malformation or deterioration of the material.	●Do not mark on the product with a pen or leave printed materials in contact with its surface. Ink marks may be absorbed and become irremovable.

Reference


Reference

Aqua Slab Measurement Result
(BKG Corrected)



Aqua Slab, a water-equivalent material, exhibits high equivalence to water across a wide energy range.

Spectral HU Curve



Aqua Slab rods containing varying iodine concentrations simulate the X-ray attenuation of contrast-enhanced blood vessels.

⚠ Caution **Don’ t mark on the phantom with pen or leave printed materials contacted on its surface. Ink marks on the phantom will be irremovable.**

- If you have any questions, please contact the dealer you purchased the product or Kyoto Kagaku Co., Ltd. below.

KYOTO KAGAKU CO.,LTD

URL: http://www.kyotokagaku.com e-mail: rw-kyoto@kyotokagaku.co.jp

Kyoto Kagaku Head Office and Factories
15 Kitanekoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN
Tel: +81-75-605-2510 Fax: +81-75-605-2519

“Worldwide Inquiries and Ordering”
Asia and other regions
Kyoto Kagaku / Overseas Sales Division
Second fl. NREG Hongo-3-chome bldg., 3-26-6 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN

All American regions
KYOTOKAGAKU AMERICA INC.
3109 Lomita Boulevard, Torrance, CA 90505-5108, USA
Tel: 1-310-325-8860 Fax: 1-310-325-8867

Europe, Russia & Africa
Kyoto Kagaku Europe GmbH
De-Saint-Exupery-Str.10, 60549 Frankfurt, Germany
Tel:+49-69-5060-28160

2026.10