
Penile Ultrasound Examination Phantom

For practical, hands-on training in penile ultrasound examination





Item No.	Code
US-26	41961-000

FEATURES

| 01 Anatomically accurate structures provide clinically realistic ultrasound images.
| 02 The flexible penis allows natural handling and positioning for scanning at various angles.
| 03 Removable phantom part allows easy cleaning of residual gel after training.

SKILLS

Ultrasound examination of the penis and penile urethra

[Recommended Devices]
Linear transducer

CASES / PATHOLOGIES / ANATOMY

Corpora cavernosa of the penis
Corpus spongiosum
Urinary tract
Blood vessels

KEY FEATURES

Anatomical details supports convincing hands-on training

DESCRIPTIONS

SPECIFICATIONS

Main unit: adult male pelvic model
Size: W34 x D33 x H24 cm (W13.4 x D13 x H9.4 in)
Material: Soft resin/Latex free

SET INCLUDES

1 Pelvic body
1 Normal unit
1 Tumor unit
1 Talcum powder
1 Storage case